



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

Quick-RNA™ Midiprep Kit

Catalog No. **R1056**

Highlights

- 10 minute method for isolating RNA (up to 1 mg) from a wide range of cell types and tissue samples.
- Clean-Spin™ column technology allows ultra clean RNA to be eluted into minimal volumes ($\geq 200 \mu\text{l}$).
- Omits the use of organic denaturants and proteases.

Contents

Product Contents.....	1
Specifications.....	1
Product Description.....	2
Buffer Preparation.....	3
Protocol.....	3
Ordering Information.....	4

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682.

Notes:

¹ Add 24 ml 100% ethanol (26 ml of 95% ethanol) to the 6 ml **RNA Wash Buffer** concentrate before use.

² For purification of RNA from whole blood, see the **Quick-RNA™ Whole Blood** (R1201).

Product Contents

Quick-RNA™ Midiprep Kit (Kit Size)	R1056 (25 Preps.)
ZR RNA Buffer	100 ml
RNA Pre-Wash Buffer	12 ml
RNA Wash Buffer¹ (concentrate)	6 ml
DNase/RNase-Free Water	10 ml
Zymo-Spin™ V-E Columns w/ Reservoir	25
Collection Tubes	50
Instruction Manual	1

Note – Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Storage Temperature – Store all components (i.e., buffers, columns) at room temperature.

Specifications

- **Sample Sources** – Cells from culture or tissue samples, buffy coat, plasma, serum, and other biological liquids. Not compatible with whole blood.²
- **Sample Size** – ~10³ to 10⁸ mammalian (animal) cells.
- **Format** – Spin column.
- **RNA Purity** – High quality RNA ($A_{260}/A_{280} >1.8$, $A_{260}/A_{230} >1.8$) suitable for all downstream RNA-based manipulations.
- **RNA Recovery** – Up to 1 mg RNA can be eluted into ≥200 µl RNase-free water allowing for a highly concentrated sample.
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored at ≤-70 °C. The addition of RNase inhibitors is highly recommended for prolonged storage.
- **Equipment Needed** – Vacuum manifold, microcentrifuge.


Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

ZYMO RESEARCH CORP.

Product Description

The **Quick-RNA™ Midiprep Kit** is an innovative product designed for the easy, reliable, and rapid isolation of up to 1 mg total RNA from cultured cells or tissue samples. The procedure combines a unique, single-step RNA extraction/binding buffer with Clean-Spin™ column technology to yield high quality RNA in about 10 minutes. The **Quick-RNA™ Midiprep Kit** allows for the efficient recovery of total RNA from 10^3 to 10^8 cells or tissue.

The method is easy: simply add the provided **ZR RNA Buffer** to extract total RNA from the cells of interest then purify the RNA using the provided **Zymo-Spin™ V-E Columns**. The result is highly-concentrated, purified RNA that is suitable for subsequent RNA-based methods including RT-PCR, hybridization, etc.



PCR amplification of β -actin transcript (353 bp fragment shown) following RNA isolation from human epithelial cells (HCT 116) with the **Quick-RNA™ Midiprep Kit**: RT-PCR (lane 1, 2), PCR negative control (RNA template; lane 3, 4).

Quick-RNA™ - Fast-Spin column and plate technology overview



	Quick-RNA™ Microprep Kit	Quick-RNA™ Miniprep Kit	Quick-RNA™ Midiprep Kit	Quick-RNA™ 96 Kit
<i>Format</i>	Zymo-Spin™ IC Column	Zymo-Spin™ IIICG Column	Zymo-Spin™ V-E Column	Silicon-A™ Plate
<i>Catalog</i>	R1050/R1051	R1054/R1055	R1056	R1052/R1053
<i>Preps</i>	50/200	50/200	25	2x 96/4x 96
<i>Input</i>	~1- 10^6 cells	~ 10^2 - 10^7 cells	~ 10^3 - 10^8 cells	~1- 10^6 cells
<i>Binding</i>	~10 μ g	~100 μ g	~1 mg	~10 μ g/well

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

Make sure guidelines are followed to ensure the RNA isolation procedure is performed in an RNase-free environment.

Clean-Spin™ column technology efficiently removes the majority of DNA during RNA purification and is satisfactory for most RNA-based applications. However, if necessary, complete removal of DNA can be achieved by performing a DNase I digestion. See the **RNA Clean & Concentrator™**.

Notes:

¹ In order to lyse samples completely, the amount of the **ZR RNA Buffer** should be adjusted (*i.e.*, more buffer can be added).

² Total RNA including small RNAs ≥ 17 nt is recovered. **Zymo-Spin™** columns may be reloaded.

³ Total RNA without small RNAs (*e.g.*, tRNAs and microRNAs) is recovered

⁴ Set vacuum source at ≥ 500 mm Hg.

⁵ To maximize RNA yield, increase the elution volume (≥ 200 μ l) and/or repeat the elution.

Buffer Preparation

Before starting, add 24 ml 100% ethanol (26 ml 95% ethanol) to the 6 ml **RNA Wash Buffer** concentrate.

Protocol

All centrifugation steps should be performed at 10,000–16,000 x g, unless specified otherwise.

1. **Sample Homogenization and Cell Lysis:** Follow an applicable method (A-D).
 - A. **Adherent Cells:** Cells can be lysed directly in the culture container by removing liquid medium and adding **ZR RNA Buffer**¹ directly to the monolayer (*e.g.*, 3 ml for 5×10^6 cells). Remove cells from culture surface by pipetting, scraping, *etc.*
 - B. **Cells in Suspension:** Pellet cells by gentle centrifugation (*e.g.*, 5 minutes at 500 x g). Remove the supernatant completely and resuspend the cell pellet in 3 ml **ZR RNA Buffer**¹. Vortex briefly.
 - C. **Solid Tissue Samples:** Add 3 ml ZR RNA Buffer¹ to fresh or frozen tissue (up to ~100 mg) and homogenize the sample (*e.g.*, using a Dounce or similar homogenizer). **Important!** To eliminate tissue debris, if present, sample lysates may be spun at ≤ 500 x g for 1 minute.
 - D. **Liquid Samples:** Add 3 volumes of **ZR RNA Buffer**¹ to every volume of sample (*e.g.*, 3 ml of buffer to 1 ml sample).

For total RNA including small RNAs ²	For total RNA without small RNAs ³
Add 1 volume ethanol (95-100%) to the lysate and mix thoroughly by pipetting.	Go directly to Step 2.

2. Transfer lysate into the **Zymo-Spin™ V-E Column with Reservoir** mounted on a vacuum manifold and start vacuum⁴.
3. Remove the reservoir and transfer column into a **Collection Tube**. Centrifuge the column for 30 seconds.
4. Add 400 μ l **RNA Pre-Wash Buffer** to the column and centrifuge for 1 minute. Discard the flow-through.
5. Add 400 μ l **RNA Wash Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through. Repeat this step.
6. Centrifuge the column for 2 minutes in the emptied collection tube to ensure complete removal of the wash buffer. Then, carefully transfer the column into an RNase-free tube (not provided).
7. Add ≥ 200 μ l **DNase/RNase-Free Water**⁵ directly to the column matrix and let stand at room temperature for 1 minute. Centrifuge for 1 minute.

The eluted RNA can be used immediately or stored at ≤ -70 °C.

ZYMO RESEARCH CORP.

Ordering Information

Product Description	Input	Binding	Catalog No.	Kit Size
Quick-RNA™ Microprep Kit	~1-10 ⁶ cells	~10 µg	R1050	50 Preps.
			R1051	200 Preps.
Quick-RNA™ Miniprep Kit	~10 ² -10 ⁷ cells	~100 µg	R1054	50 Preps.
			R1055	200 Preps.
Quick-RNA™ Miniprep Plus Kit	~10 ² -10 ⁷ cells	~100 µg	R1057T	10 Preps.
			R1057	50 Preps.
			R1058	200 Preps.
Quick-RNA™ Midiprep Kit	~10 ⁶ -10 ⁸ cells	~1 mg	R1056	25 Preps.
Quick-RNA™ 96 Kit	~1-10 ⁶ cells	~10 µg/well	R1052	2x 96 Preps.
			R1053	4x 96 Preps.

For Individual Sale	Catalog No.	Amount
ZR RNA Buffer	R1020-1-50	50 ml
	R1020-1-100	100 ml
	R1020-1-200	200 ml
RNA Pre-Wash Buffer	R1020-2-12	12 ml
	R1020-2-25	25 ml
	R1020-2-50	50 ml
	R1020-2-100	100 ml
RNA Wash Buffer (concentrate)	R1003-3-6	6 ml
	R1003-3-12	12 ml
	R1003-3-24	24 ml
	R1003-3-48	48 ml
Zymo-Spin™ V-E Columns w/ Reservoir	C1029-25	25
Collection Tubes	C1001-50	50
	C1001-500	500
	C1001-1000	1000
DNase/RNase-Free Water	W1001-1	1 ml
	W1001-4	4 ml
	W1001-6	6 ml
	W1001-10	10 ml

ZYMO RESEARCH CORP.

RNA MADE SIMPLE



ZYMO RESEARCH
The Beauty of Science is to Make Things Simple

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com