

Easy DNA Isolation from FFPE Samples


Quick-DNA™ FFPE Kit

- **Quick & Easy:** Rapid dewaxing procedure (no xylene necessary).
- **Ultra-Pure:** Ready for qPCR, Next-Gen Sequencing, arrays, etc.
- **Highest Yield:** Recover 6x more DNA.

Simple Workflow




The Highest Recovery



DNA isolated using the Quick-DNA™ FFPE Kit consistently yielded more usable DNA, indicated by lower Ct values, than the Supplier Q kit. Equivalent amounts of DNA isolated using the Zymo Research and Supplier Q procedures were used for real time PCR analysis.

Size Selection Built In



The Quick-DNA™ FFPE Kit selectively isolates DNA >50 bp or >500 bp. Equivalent amounts of DNA were resolved on a 1% agarose gel stained with EtBr.

Product	Cat. No.	Size	Binding Capacity	Minimum Elution	Input
Quick-DNA™ FFPE Kit	D3067	50 preps	25 µg	25 µl	≤ 25 mg