


# Easy & Fast Protein Purification from Crude Lysates

## MBP-Spin™ Protein Miniprep Kit

- **Fast & Simple:** Purify MBP-tagged proteins from cell-free extracts using a spin-column in  $\geq 6$  minutes.
- **High-Quality:** Prepare pure proteins for small-scale studies.
- **Convenient:** No special instrumentation needed other than a microcentrifuge.



The MBP-Spin™ Protein Miniprep Kit purifies high-quality MBP-tagged proteins directly from a spin-column. N-terminal MBP-tag fusion proteins were expressed in *E. coli* cells, and the cell extracts as well as the proteins purified using the MBP-Spin™ Protein Miniprep Kit were analyzed by SDS-PAGE on a 4-20% gel and stained with InstantBlue™; (MBP-Calmodulin 55 kDa, MBP-BDV-Phosphoprotein 65 kDa, MBP-GFP 69 kDa).

## Simple & Fast



Product	Cat. No.	Size
MBP-Spin™ Protein Miniprep Kit	P2006 P2007	10 preps 50 preps