

PCR Machine Manufacturer	Model Number	Compatible with SKU
ABI Applied Biosystems	Quant Studio 3	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02
ABI Applied Biosystems	Quant Studio 5	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02, PCR-384-ABI
ABI Applied Biosystems	Quant Studio 6	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02, PCR-384-ABI
ABI Applied Biosystems	Quant Studio 7	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02, PCR-384-ABI
ABI Applied Biosystems	Quant Studio 12K	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02, PCR-384-ABI
ABI Applied Biosystems	7000	PCR-HS-FL02,PCR-HS-RR02 ,PCR-HS-02
ABI Applied Biosystems	7300	PCR-HS-FL02,PCR-HS- RR02,PCR-HS-02
ABI Applied Biosystems	7500	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02
ABI Applied Biosystems	7900HT	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS- 02,PCR-384-ABI
ABI Applied Biosystems	9700	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-02,PCR-384-ABI
ABI Applied Biosystems	9800	PCR-HS-RR01
ABI Applied Biosystems	2700	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-02
ABI Applied Biosystems	2720	PCR-HS-FL02, PCR-HS- RR02,PCR-HS-02
ABI Applied Biosystems	StepOne Plus	PCR-HS-RR01
ABI Applied Biosystems	Automated Thermal Cycler	PCR-HS-RR02 PCR-384-ABI
ABI Applied Biosystems	MiniAmp	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-02
ABI Applied Biosystems	MiniAmp Plus	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-02
ABI Applied Biosystems	SimpliAmp	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-02
ABI Applied Biosystems	Veriti	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02, PCR-384-ABI
ABI Applied Biosystems	VeritiDX	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02, PCR-384-ABI
ABI Applied Biosystems	ViiA7	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-RR01,PCR-HS-02,

**MOLECULAR BIOLOGY PRODUCTS, INC**

Phone: 416.549.8100 | [www.mbpinc.net](http://www.mbpinc.net) | [orders@mbpinc.net](mailto:orders@mbpinc.net)

PCR Machine Manufacturer	Model Number	Compatible with SKU
		PCR-384-ABI
ABI Applied Biosystems	ProFlex	PCR-HS-FL02,PCR-HS-RR02 PCR-HS-02
Agilent	AriaMX	PCR-NS-02,PCR-FS-02
Agilent	Gradient Cyler	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Agilent	Robocycler	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Agilent	Mx3000P	PCR-NS-02,PCR-HS-02
Agilent	Mx3005P	PCR-NS-02,PCR-HS-02
Agilent	Mx4000	PCR-NS-02,PCR-HS-02
Agilent	SureCycler8800	PCR-NS-02,PCR-NS-01,PCR-384-ABI
Analytik Jena/Biometra	FlexCyler2	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Analytik Jena/Biometra	qTower	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Analytik Jena/Biometra	SpeedCyler2	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Analytik Jena/Biometra	T1	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Analytik Jena/Biometra	Tadvanced	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Analytik Jena/Biometra	TGradient	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Analytik Jena/Biometra	TOne	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Analytik Jena/Biometra	TProfessional	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Analytik Jena/Biometra	TRobot	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Axygen	MaxyGene II	PCR-HS-02
Bio-Rad	C1000	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Bio-Rad	C1000 Touch	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Bio-Rad	CFX Connect	PCR-NS-01,PCR-FS-02,PCR-384-ABI
Bio-Rad	CFX Touch	PCR-NS-01,PCR-FS02,PCR-

**MOLECULAR BIOLOGY PRODUCTS, INC**

Phone: 416.549.8100 | [www.mbpinc.net](http://www.mbpinc.net) | [orders@mbpinc.net](mailto:orders@mbpinc.net)

PCR Machine Manufacturer	Model Number	Compatible with SKU
		384-ABI
Bio-Rad	Chromo4	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Bio-Rad	DNA Engine	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Bio-Rad	Opticon	PCR-NS-01,PCR-FS-02
Bio-Rad	iCycler	PCR-NS-02,PCR-HS-02
Bio-Rad	iCycler 1Q	PCR-NS-02,PCR-HS-02
Bio-Rad	iQ5	PCR-NS-02,PCR-HS-02
Bio-Rad	PTC-200	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Bio-Rad	S-1000	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Bio-Rad	T100	PCR-NS-02,PCR-HS-02
Bioer/Bulldog	Gene Touch	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Bioer/Bulldog	Gene Touch Plus	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Bioer/Bulldog	LifeECO	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Bioer/Bulldog	LifePro	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Bioer/Bulldog	LifeTouch	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Bruker/Hain Lifesciences	FluoroCycler XT	PCR-FS-01W
Bruker/Hain Lifesciences	GTQ-Cyler 96	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Cole Parmer/Techne	Flexigene	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Cole Parmer/Techne	Genius	PCR-NS-02,PCR-NS-01,PCR-HS-02
Cole Parmer/Techne	Prime	PCR-NS-02,PCR-NS-01,PCR-HS-02
Cole Parmer/Techne	Prime Elite	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Cole Parmer/Techne	Prime Q	PCR-NS-01,PCR-FS-02
Cole Parmer/Techne	Quantica	PCR-NS-01, PCR-FS-02
Cole Parmer/Techne	TC-4000	PCR-NS-02, PCR-NS-01,PCR-HS-02,PCR384-ABI
Cole Parmer/Techne	TC-412	PCR-NS-02,PCR-NS-

**MOLECULAR BIOLOGY PRODUCTS, INC**

Phone: 416.549.8100 | [www.mbpinc.net](http://www.mbpinc.net) | [orders@mbpinc.net](mailto:orders@mbpinc.net)

PCR Machine Manufacturer	Model Number	Compatible with SKU
		01,PCR-HS-02,PCR-384-ABI
Cole Parmer/Techne	TC-5000	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-384-ABI
Cole Parmer/Techne	TC-512	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Cole Parmer/Techne	TC-PLUS	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
DLAB	TC1000-G	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Eppendorf	Master ep Gradient	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Eppendorf	Master ep realplex	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Eppendorf	Mastercycler Gradient	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Eppendorf	Mastercycler Nexus	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Eppendorf	Mastercycler Pro	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Eppendorf	Mastercycler X50	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Eastwin	ETC811	PCR-NS-02,PCR-HS-02
Eastwin	ETC821	PCR-NS-02,PCR-HS-02
Esco	Aeris	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Esco	Swift MaxPro	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Genestorm/Agilegene	GS1	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Genestorm/Agilegene	GS2	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Genestorm/Agilegene	GS4	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-

**MOLECULAR BIOLOGY PRODUCTS, INC**

Phone: 416.549.8100 | [www.mbpinc.net](http://www.mbpinc.net) | [orders@mbpinc.net](mailto:orders@mbpinc.net)

PCR Machine Manufacturer	Model Number	Compatible with SKU
		02,PCR-384-ABI
Kyratec	SuperCycler Trinity	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
Labnet	MultiGene	PCR-NS-02,PCR-HS-02
Labnet	OptiMax	PCR-NS-02,PCR-HS-02
MWG Biotech	Primus	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
MWG Biotech	The Q LifeCycler	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
PeqLab/VWR	PeqSTAR2X	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02
PeqLab/VWR	peqSTAR X	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
PeqLab/VWR	Uno	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Roche	LightCycler 480	PCR-HS-01W,PCR-384-ROC
Roche	LightCycler 96	PCR-HS-01W
SensoQuest	Labcyler	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
ThermoFisher	Arktik	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
ThermoFisher	Hybaid MBS	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
ThermoFisher	Hybaid Px2	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
ThermoFisher	Hybaid PxE	PCR-NS-02,PCR-NS-01,PCR-HS-02,PCR-FS-02,PCR-384-ABI
Takara	TP3000	PCR-NS-02,PCR-NS-01,PCR-HS-02